

BONDED PHASES FOR HPLC AND THEIR ABBREVIATIONS

HPLC 固定相簡介

Phase	Description	Chemical Structure	Phase	Description	Chemical Structure
Si	Silica Classic normal phase material. Suitable for separating polar non-ionic organic compounds.	$\begin{array}{c} \\ -\text{Si}-\text{OH} \\ \end{array}$	CN	CPS, PCN, Cyano, Cyanopropyl, Nitrile Can be employed as either a reversed phase or normal phase material. Slightly polar, unique selectivity for polar compounds in both reversed phase and normal phase modes. Equilibrates very rapidly, suitable for gradient separations. Useful for many pharmaceutical applications (e.g. tricyclic antidepressants).	$\begin{array}{c} \\ -\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{CN} \\ \end{array}$
C1	TMS, SAS, Trimethyl Reversed phase material. Unique selectivity for polar and multifunctional compounds. Least retentive of all alkyl group bonded phases for non-solutes.	$\begin{array}{c} \\ -\text{Si}-\text{CH}_3 \\ \end{array}$	C2	RP-2, Dimethyl Reversed phase material. Less retentive than C4, C8, or C18. More retentive than C1.	$\begin{array}{c} \\ -\text{Si}-\text{C}_2\text{H}_5 \\ \end{array}$
C2	RP-2, Dimethyl Reversed phase material. Less retentive than C4, C8, or C18. More retentive than C1.	$\begin{array}{c} \\ -\text{Si}-\text{C}_2\text{H}_5 \\ \end{array}$	NH ₂	APS, Amino, Amino Propyl Silyl Can be employed as reversed phase, normal phase, or weak anion exchange material. Reversed phase: useful for separating carbohydrates. Normal phase: alternative selectivity to silica, not deactivated by small amounts of water. Ion exchange: weak anion-exchanger when used with buffers. Separates anions and organic acids.	$\begin{array}{c} \\ -\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{NH}_2 \\ \end{array}$
C3	Propyl Reversed phase material. Used in Hydrophobic Interaction Chromatography (HIC) of proteins and peptides.	$\begin{array}{c} \\ -\text{SiC}_3\text{H}_7 \\ \end{array}$	NO ₂	Nitro Normal phase material. Separates aromatic compounds and compounds with double bonds.	$\begin{array}{c} \\ -\text{Si}-\text{NO}_2 \\ \end{array}$
C4	Butyl Reversed phase material. Useful for HIC and ion-pairing chromatography. Offers less retention than C8 and C18 phases for non-polar solutes. When bonded to 300Å silica, it is an ideal phase for analysing large proteins and hydrophobic peptides.	$\begin{array}{c} \\ -\text{Si}-\text{C}_4\text{H}_9 \\ \end{array}$	OH	Diol, Glycerol Can be employed as either a reversed phase or normal phase material. Reversed phase: used for Gel Filtration Chromatography (GFC) of proteins and peptides. Normal phase: similar selectivity to silica, not deactivated by small amounts of water.	$\begin{array}{c} \\ -\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{OCH}_2\text{CH}(\text{OH})\text{CH}_2\text{OH} \\ \end{array}$
C5	Pentyl Reversed phase material. Useful for reversed phase separation of hydrophobic proteins and peptides when bonded to 300Å silica. More hydrolytically stable than C4.	$\begin{array}{c} \\ -\text{SiC}_5\text{H}_{11} \\ \end{array}$	SAX	SB, Quaternary amine, Strong Base Ion-exchange material. Strong anion-exchangers (basic) are useful for separating nucleotides, nucleosides, and organic acids.	$\begin{array}{c} \\ -\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{N}^+(\text{CH}_3)_3 \\ \end{array}$
C6	Hexyl Reversed phase material. Useful for ion-pairing chromatography. Less retentive than C8 and C18 phases.	$\begin{array}{c} \\ -\text{Si}-\text{C}_6\text{H}_{13} \\ \end{array}$	SCX	SA, Sulfonic Acid, Strong Acid Ion-exchange material. Strong cation-exchangers (acidic) are useful for separating organic bases.	$\begin{array}{c} \\ -\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2\text{SO}_2\text{OH} \\ \end{array}$
C8	MOS, RP-8, LC8, Octyl Reversed phase material. Similar selectivity to C18 but less retentive. Wide applicability (e.g., pharmaceuticals, nucleosides, steroids). When bonded to 300Å silica, it is an ideal phase for peptides, peptide mapping, and small hydrophilic proteins.	$\begin{array}{c} \\ -\text{Si}-\text{C}_8\text{H}_{17} \\ \end{array}$	WAX	DEAE, Diethylaminoethyl, Weak Base PEI, Polyethyleneimine Ion-exchange material. Weak anion-exchangers (acidic) are most useful for analyzing acidic proteins and peptides.	$\begin{array}{c} \\ -\text{Si}-\text{CH}_2\text{CH}_2\text{N}(\text{CH}_2\text{CH}_3)_2 \\ \\ \\ -\text{Si}-\text{CH}_2\text{CH}_2\text{NH} \\ \end{array}$
C18	ODS, RP-18, LC18, Octadecyl Classic reversed phase material. Most retentive for nonpolar solutes. Excellent for ion-pairing chromatography. Wide applicability (e.g. nucleosides, nucleotides, steroids, pharmaceuticals, vitamins, fatty acids, environmental compounds). When bonded to 300Å silica, this phase is perfect for separating small hydrophilic peptides.	$\begin{array}{c} \\ -\text{Si}-\text{C}_{18}\text{H}_{37} \\ \end{array}$	WCX	CM, Carboxymethyl, Weak Acid Ion-exchange material. Weak cation-exchangers (basic) are most useful for analyzing basic proteins and peptides.	$\begin{array}{c} \\ -\text{Si}-\text{CH}_2\text{COOH} \\ \end{array}$
Phenyl	C6H5 Reversed phase material. Unique selectivity. Useful for analyzing aromatic compounds. When bonded to 300Å silica, this phase is useful for HIC.	$\begin{array}{c} \\ -\text{Si}-\text{CH}_2\text{CH}_2\text{CH}_2-\text{C}_6\text{H}_5 \\ \end{array}$			