
Based on Stepwise Pressure Programming Method

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FEATURES

- STEPWISE PRESSURE PROGRAMME
Pressure Range, 0~730 kg/cm²
- HIGH FLOW, PULSELESS
300ml/min (Unloaded)
- RESISTANCE TO ALL SOLVENTS
Solvent-Switching mode
- UPWARD PACKING
Modified U Slurry Reservoir
- UNIQUE PACKING METHOD
With use of low viscous, self-dispersible
Chemco Slurry Mix

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
 **Chemco**

HERE IS YOUR ANSWER TO PACK HIGH EFFICIENCY HPLC COLUMNS WITH OUR SPECIFIC KNOWHOWS, WHICH WILL BE PRECISELY TRANSFERRED TO THE USERS OF OUR APPARATUS

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TEL 06 (6375) 0395 (代表) FAX 06 (6372) 0533

TEL 03 (3444) 4201 (代表) FAX 03 (3444) 4238

CHEMCO AUTO PRESSURE PROGRAMMER MODEL 124PP

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■ NO HUMAN ERRORS FOR COLUMN-PACKING (カラム充填の個人差排除)

Column Packing Pressure can be automatically programmed with pre-set profile based on our complete knowhows for available silica-based packing media. Our conventional manual step programming rates are easily converted to auto programming ones; maximum pressure needed V.S. time consumed.

Besides, this new programmer also controls time hold at pre-set maximum pressure, and immediately after this hold time ends, the pressure drops to atmosphere during 50 seconds.



Then, after for a while, you can dejoint the packed tube from the slurry reservoir, and set the end-reducer for HPLC columns. That's all you do for column packing; first place the packing slurry into the slurry-reservoir, set the maximum pressure, program rate, hold time, and then you only push the start switch. All is in your simple operation for HPLC column-packing with this new device, M-124PP.

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
Based on Automatic Pressure Programming

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FEATURES

We have placed, on the market, Model 124PP Auto-pressure Programmer as an optional accessory for our conventional Model 124A. This unique accessory enables HPLC-column packing unattended, and constant in column performance. The packing procedure is very simple that: ① Place the packing slurry into the slurry reservoir, ② Set the programme rate and the maximum pressure needed, and ③ Push the start switch; that's all you have to do!

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CO., LTD.
Chemco Scientific Co., Ltd.