# chromatography products





TRANSGENOMIC<sup>®</sup> BIOCONSUMABLES<sup>™</sup>



Transgenomic is a global company focused on providing you the best separations technology with the highest reproducibility possible. We understand the quality of your results depends on us. Our entire team is dedicated to supporting you in your scientific quest.

The separations products we provide are based on our many years of experience in developing and manufacturing polymer chemistries for liquid chromatography. This vast experience and knowledge base continues to help us build on our strong tradition of providing the best products for your research and quality control needs.

Collin C.S.

Collin D'Silva Chief Executive Officer

## TABLE OF Contents

Application Guide	2
Amino Acid Analysis	3
Carbohydrate Analysis	7
Organic Acid Analysis	21
Polymeric Reversed Phase	32
Ion Chromatography	38
Guard-Disc <sup>®</sup> Protection System	42
Solid Phase Extraction	44
Bulk Polymeric Resin	46
Buffers	47
Hardware	48
Index Ir	side backcover





## APPLICATION SELECTION Guide

Amino Acids	Protein Hydrolysates	AMINOSep AA511 AMINOSep AA911 Na⁺ Column for 63/7300 Systems Na⁺ Column for System Gold
	Physiological Fluids	Li⁺ Column for 63/7300
Carbohydrates	Monosaccharides Disaccharides Sugar Alcohols	CARBOSep CHO-620 CARBOSep CHO-682 CARBOSep CHO-820 CARBOSep CHO-6110H CARBOSep USP L-19 CARBOSep COREGEL-87C CARBOSep COREGEL-87P CARBOSep COREGEL-87H CARBOSep COREGEL-87MM
	Oligosaccharides, Corn Syrup,, Sugar Polymers	CARBOSep COREGEL-42Ag CARBOSep CHO411 CARBOSep CHO611 CARBOSep COREGEL-87K CARBOSep COREGEL-87N
Organic Acids	Sugar Alcohols Organic Acids	ICSep Ion-300 ICSep COREGEL-87H ICSep COREGEL-107H ICSep ORH-801 ICSep WA-1 Wine Analysis Column ICSep Ion-310 ICSep ARH-601 ICSep COREGEL 64H
Proteins/Peptides	Reversed Phase	RPSep ACT-1 C18 RPSep PRX-1 RPSep PolyRP C0
DNA, RNA, Oligonucleotides	Reversed Phase	RPSep PRX-1

# AMINO ACID Analysis

#### Transgenomic Columns for Amino Acid Analysis

Ion-exchange chromatography is a popular technique for the analysis of amino acids because both retention times and quantification are highly reproducible regardless of the sample matrix. This unique matrix insensitivity is important when comparing results from different patients or batches of protein hydrolysate.

Amino acids are zwitterions; at low pH, they are positively-charged and are bound to the resin by their attraction to the negativelycharged ion-exchange sites. Almost all the contaminants, i.e. matrix, are eluted at the void. The amino acids are then selectively eluted by increasing the pH and salt concentration with different buffers. With few exceptions, the order of elution follows the isoelectric point of the amino acids, i.e. acidic amino acids first, then neutral and basic. Because the separation and the ensuing post-column reaction of amino acids are devoid of contaminants, amino acid analyses via ion-exchange chromatography are highly reproducible.

#### Features

The key features of the Transgenomic cation-exchange columns are:

- Polymeric Substrate
- High efficiency
- High resolution
- Reproducibility lot-to-lot and column-to-column
- Rugged
- Available for both physiological and protein hydrolysate amino acids

Amino acid columns are subjected to many different types of samples (blood, urine, growth media, animal feed, wine, etc.) and often they are introduced with minimum sample preparation. Therefore this variety of matrix challenges all but the most rugged ion-exchange columns. Transgenomic columns use polystyrene/ divinylbenzene copolymers and are stable in the pH range of 0 to 14; they are temperature stable and very rugged. The Transgenomic amino acid columns have been shown to last for thousands of runs without cleaning. Because Transgenomic manufacture the polymers and pack the columns, lot-to-lot and column-to-column reproducibility is excellent (retention times vary by less than 1%). Available for both routine hydrolysate analysis as well as complex physiological fluids, Transgenomic amino acid columns have been designed to provide the highest efficiency and highest resolution of any ion-exchange amino acid columns on the market.



#### **Oxidized Hydrolysate Standards**

#### **Analysis Conditions:** Sample: Column: Transgenomic Sodium Column for 6300 1.L-Cysteic Acid Flow rate: 0.233 mL/min 2. Methionine Sulfoxide Temperature: 48-70-77°C 3. L-Aspartic Acid Pressure: 655 PSIG 4. Methionine Sulfone Detection: Fluorescence 19 20 6 8 1011 13 Injection: 20 µL 23 12 17 3 5 22 16 15 2 21 30 50 60 ò 10 20 40

5. L-Threonine 6.L-Serine 7.L-Glutamic Acid 8. Glycine 9. L-Alanine 10. L-Valine 11.L-Methionine 12. L-Isoleucine 13.L-Leucine 14. Norlfuline 15. L-Tyrosine 16. L-Phenylalanine 17. Glucosamine 18. Galactosamine 19. L-Histidine 20. L-Lysine 21.Tryptophan 22. Ammonia 23. L-Arginine

#### Physiological Fluid Amino Acids

#### **Analysis Conditions:** 20.Valine Sample: Column: Transgenomic Lithium Column for 6300 21. Cystine 1. Phosphoserine Flow rate: 0.333 mL/min 22. Methionine 2. Taurine Temperature: 32.5-63-80°C 23. Cystathionine 3. Phosphoethanolamine Pressure: 1200 PSIG 24. Isoleucine 4. Urea Detection: UV 25. Leucine 5. Glucosaminic Acid Injection: 20 µL 26. Tyrosine 6. Aspartic Acid 27. Phenylalanine 7. Hydroxyproline 28. $\beta$ -Alanine 8.Threonine 29. $\beta$ -Aminoisobutyric Acid 9. Serine . 30. Homocystine 10. Asparagine 31. $\gamma$ -Aminobutyric Acid 11. Glutamic Acid 21 32. Ethanolamine 12. Glutamine 13. Sarcosine 33. Ammonia 14. $\alpha$ -Aminoadipic Acid 34. Hydroxylysine 35. allo-Hydroxylsine 15. Proline 36. Aminoethylcysteine 16. Glycine 37. Ornithine 17. Alanine 2 23 6 38. Tryptophan 18. Citrulline 39. Lysine 19. $\alpha$ -Amino- $\eta$ -Butyric Acid 40.1-Methylhistidine 30 41. Histidine 22 42.3-Methylhistidine 5 43. Anserine 3 362-44. Carnosine 45. Arginine 25 41 Muh | 110 min. 10 20 30 40 50 60


#### **Amino Acid in Red Wine**

#### Analysis Conditions:

Column: Transgenomic Sodium Column for 6300 Flow rate: 0.233 mL/min Temperature: 48-70-77°C Pressure: 575 PSIG Detection: Fluorescence Injection: 20 µL

#### Sample:

1. Cysteic Acid
2. ASP
3.MT02
4.THR
5.GLU
6. GLY
7. ALA
8. MET
9. Glucosamine
10. Galactosamine
11.HIS
13.LYS
14.NH 3
15. A R G




#### **Amino Acid in Urine**

#### Analysis Conditions:

Column: Transgenomic Lithium Column for 6300 Flow rate: 0.333 mL/min Temperature: 32.5-63-80°C Pressure: 1200 PSIG Detection: Fluorescence Injection: 20 µL

Sample:	16.T R P
1.PER	17.EIN
2.TAU	18.NH 3
3.PETN	19.0 R N
4.THR	20. L Y S
5.GLU	21.1 ME-HIS
6. GLY	22.HIS
7.ALA	23.3 ME-HIS
8. Met	24. A N S
9.CYST	25.CARN
10.ILE	26. A R G
11.LEU	
12.TYR	
13.PHE	
14. BALA	
15.BABA	



#### Transgenomic Lithium Amino Acid Column

#### (4 x 100 mm)

#### P/N AAA-99-6311

- Designed for use with the Beckman Coulter® 6300 and 7300 Amino Acid Analyzers using either the Beckman or Pickering Lithium buffer systems
- The Lithium column is ideal for Physiological amino acid analysis
- Highly efficient 6 micron particle size

#### **AMINOSep Lithium Guard Kit**

P/N AAA-99-2311

#### AMINOSep Lithium Guard Cartridge – 2/PK

P/N AAA-99-1311

#### Transgenomic Sodium Amino Acid Column

#### (4 x 120 mm) P/N AAA-99-6312

- P/N AAA-99-6312
- Designed for use with the Beckman Coulter 6300 and 7300 Amino Acid Analyzers using either the Beckman Coulter or Pickering Sodium buffer systems
- The Sodium column is ideally suited for routine hydrolysate analysis
- Extremely rugged polymer

#### **AMINOSep Sodium Guard Kit**

P/N AAA-99-2312

#### AMINOSep Sodium Guard Cartridge – 2/PK P/N AAA-99-1312

#### Transgenomic Sodium Sodium Amino Acid Column for Use with System Gold

#### (4 x 200 mm)

#### P/N AAA-99-6310


- Designed for use with the Beckman Coulter System Gold Amino Acid Analyzer
- This Sodium cation exchange column is ideal for the separation of hydrolysate amino acids.


#### **AMINOSep Sodium Guard Kit**


#### P/N AAA-99-2312


### AMINOSep Sodium Guard Cartridge – 2/PK

P/N AAA-99-1312









#### AMINOSep AA-511 Sodium Column


(4.6 x 150mm) P/N AAA-99-7554

#### AMINOSep GC-511 Guard Kit

P/N AAA-99-2354

#### AMINOSep GC-511 Guard Cartridge – 2/PK

P/N AAA-99-1354



#### AMINOSep AA-511 High Speed Sodium Column


(4.6 x 120mm) P/N AAA-99-6554

#### AMINOSep GC-511 Guard Kit

P/N AAA-99-2354

#### AMINOSep GC-511 Guard Cartridge – 2/PK

P/N AAA-99-1354



## CARBOHYDRATE Analysis

#### **CARBOSep** Columns

Transgenomic manufactures a line of polymeric columns for carbohydrate analysis called CARBOSep columns. CARBOSep columns employ a technique called ligand-exchange chromatography for the separation of monosaccharides, disaccharides and oligosaccharides up to 15 glucose units long.

The principle behind ligand exchange is that each of the hydroxyls on a sugar molecule carry a very slight negative charge. The hydroxyl group on the anomeric carbon can be deprotonated and have a strong negative charge. It is the interaction between these negative charges on the sugar molecule and the positive charge contributed by the metal ion secured to the resin surface that causes the sugars to be retained and thus separated.

Ligand exchange resins are highly sulfonated cation exchange resins that have group 1, 2 or transition series metals loaded on. The sulfonic acid groups on the resin tightly hold the metal ions via an ionic attraction so that it is not released during analysis or through the life of the column. It is this metal ion that provides the positive charge that interacts with the negative charge on the sugar. During analysis, the carbohydrates are introduced onto the column. The sugars are attracted to the metals via an ionic interaction thus they become weakly bound to the metal ion on the resin. Water will also have a weak ionic interaction with the metals on the column, so the water will exchange with the sugars on the metal sites. This ionic adsorption and desorption occurs for the sugars through the column. Since the ionic charge is different for every sugar, separation of the sugars occurs.

Selectivity is easily controlled by resin type, metal selected, and other factors such as temperature and mobile phase. CARBOSep columns are provided in a large variety of resin types and metals to provide selectivities that meet your separation needs.

#### Selectivity Chart for Carbohydrate Columns

Compound	CHO-620 (units in minutes)	<b>CHO-611</b> (units in minutes)	CHO-682 (units in minutes)	<b>COREGEL 87H</b> (units in minutes)	COREGEL 87P (units in minutes)	<b>COREGEL 87N</b> (units in minutes)	COREGEL 87K (units in minutes)	COREGEL 87 (units in minutes
Arabinose	10.64	11.08	23.95	12.08	16.32	12.64	14.72	13.92
Digitoxose	10.26	10.18	21.95	_	15.48	11.40	12.32	14.19
Fructose	10.07	10.33	25.84	11.25	16.96	11.61	13.31	13.63
Fucose	10.57	10.96	24.16	12.80	16.44	12.34	14.39	13.82
Galactose	9.58	10.22	22.32	11.12	15.16	11.44	13.36	13.82
Glucose	8.72	9.53	19.14	10.57	13.38	10.72	12.55	11.17
Mannose	9.79	10.27	25.50	11.13	16.76	11.57	13.74	12.76
Rhamnose	9.64	9.88	22.56	11.94	15.26	11.08	12.83	12.86
Sorbose	9.50	9.93	22.38	10.08	15.24	11.08	12.66	12.86
Tagatose	11.53	10.29	-	11.15	20.80	11.36	12.82	16.46
Xylose	9.56	10.34	20.64	11.32	14.42	11.77	13.69	12.32
Cellobiose	6.65	7.17	15.58	8.43	10.98	7.90	9.26	8.94
Lactose	7.01	7.51	17.37	8.77	11.84	8.18	9.63	9.44
Lactulose	7.57	7.85	20.70	9.00	13.24	8.48	10.08	10.17
Melibiose	6.99	7.46	17.63	8.56	12.02	8.19	9.72	9.36
Trehalose	6.70	7.14	15.98	8.64	11.20	7.85	9.02	9.07
Sucrose	6.76	7.27	15.70	_	11.10	7.99	9.11	9.09
Maltose	6.89	7.37	16.61	8.57	11.54	8.08	9.48	9.17
Ribitol	10.94	10.13	30.72	12.44	20.44	11.26	11.84	15.55
Arabitol	12.32	10.52	39.82	12.65	25.24	11.64	12.10	18.36
Galactitol	13.05	10.23	52.43	11.80	31.60	11.15	11.61	20.46
Myo-inositol	10.82	11.01	35.58	11.02	20.06	12.48	14.08	14.27
Lactitol	8.55	7.87	33.23	9.26	19.50	8.45	9.34	12.17
Maltitol	8.54	7.68	30.38	9.00	17.76	8.28	9.06	12.22
Mannitol	11.84	9.90	40.03	11.66	24.98	10.81	11.42	17.81
Sorbitol	13.64	10.38	56.56	11.77	33.40	11.32	11.86	21.34
Xylitol	13.93	11.01	51.15	12.82	31.10	12.16	12.64	21.30
Amiprylose	4.50	4.20	-	6.86	9.46	5.74	6.42	7.68
Melezitose	5.78	6.01	13.85	-	13.08	6.81	7.82	8.20
Maltotriose	5.91	6.22	15.17	7.72	10.54	6.98	8.16	8.28
Raffinose	5.86	6.10	14.40	-	10.22	6.88	7.92	8.24
Stachyose	5.28	5.39	13.41	-	9.58	6.33	7.28	7.77
Maltotetrose	5.37	5.54	14.07	7.30	9.84	6.42	7.46	7.80
Maltopentose	5.00	5.08	13.08	7.10	9.34	6.11	7.02	7.53
Maltohexose	4.78	4.87	12.24	7.00	8.80	5.94	6.74	7.38
Maltoheptose	4.66	4.60	11.74	6.96	8.52	5.84	6.61	7.28
Nitrate	4.50	4.20	10.30	6.85	8.40	5.70	6.40	7.30

Mobile Phase: 100% water. • Flow rate: 0.5 mL/minute. • Temperature: 90°C

#### **Carbohydrate Columns Specifications Chart**

Column	Application	Form	Particle Size (µm)	Typical Mobile Phase	Recom'd Rate Flow (mL/min)	Recom'd Temp (°C)
CARBOSep CHO-411	oligosaccharides up to DP10, corn syrup, molasses	sodium	20	water	0.4	75
CARBOSep CHO-611	oligosaccharides up to DP5	sodium	10	water	0.5	90
CARBOSep CHO-6110H	mono and oligosaccharides w/ PAD detection	sodium	10	sodium hydroxide	0.5	90
CARBOSep CHO-620	high fructose corn syrup, mono-, di-, trisaccharides and sugar alcohols	calcium	10	water	0.5	90
CARBOSep CHO-682	mono and disaccharides, sucrose, maltose lactose	lead	7	water	0.4	80
CARBOSep CHO-820	simple sugars, sugar alcohols	calcium	8	water	0.5	90
CARBOSep COREGEL 87C	mono and disaccharides	calcium	9	water	0.6	85
ICSep COREGEL 87H1	fast analysis of organic acids, alcohols, sugar mixtures	hydrogen	9	sulfuric acid	0.6	85
ICSep COREGEL 87H3	organic acids, alcohols, sugar mixtures	hydrogen	9	sulfuric acid	0.6	85
CARBOSep COREGEL-42Ag	oligosaccharides up to DP11	silver	20	water	0.4	75
CARBOSep COREGEL 87K	beet sugar, cane sugar, corn syrup, molasses	potassium	8	water	0.6	85
CARBOSep COREGEL 87N	beet sugars, mono and oligosaccharides	sodium	8	water	0.6	85
CARBOSep COREGEL 87P	pentose, hexose, monosaccharides, alcohols	lead	8	water	0.8	85
CARBOSep USP L19	USP L-19 specifications for separation of sorbitol and mannitol	calcium	9	water	0.2	30
CARBOSep COREGEL-87MM	mono, di, and trisaccharides, and sugar alcohols	calcium/sodium	8	water	0.5	85
ICSep ION300	glucose and fructose in organic acid mixtures	hydrogen	8	sulfuric acid	0.4	70
ICSep ION310	grape must analysis	hydrogen	8	sulfuric acid	0.8	50

Mobile Phase: 100% water. • Flow rate: 0.5 mL/minute. • Temperature: 90°C


#### **Separation of Carbohydrates with PAD**

#### **Analysis Conditions:**

Column: CHO-6110H Eluent: 0.015N NaOH Flow rate: 0.6 mL/min Temperature: 85°C Detection: PAD Injection: 5 µL


#### Sample:

1. Sucrose ( 500 ppm) 2. Glucose (250 ppm) 3. Arabinose (250 ppm)



#### Separation of Carbohydrate Standards

#### Analysis Conditions:




#### Separation of Blocked Carbohydrates

#### Analysis Conditions:


Column: CHO-6110H Eluent: 0.01 N NaOH Flow rate: 0.5 mL/min Temperature: 85°C Detection: RI Injection: 10 µL

#### Sample: 1 mg/ml each,

1. Monoacetone xylofuranose 2. Diacetone xylofuranose



#### Separation of Mannitol and Sorbitol for USP-L-19



13


#### Separation of Sugars in Apple juice

#### **Analysis Conditions:**

Column: CHO-820 (7.8 mm x 300) Eluent: Distilled Water Flow rate: 0.5 mL/min Temperature: 90°C Pressure: 50 Bar Detection: RI Range 16x Injection: 20 µL

#### Sample: **Apple Juice Diluted**

- 1:9 with DI Water
- 1. Sucrose
- 2. Glucose
- 3. Fructose
- 4. Sorbitol



3


#### **Apple Juice**

#### **Analysis Conditions:**

Column: CHO-620 Eluent: H<sub>2</sub>O Flow rate: 0.5 mL/min Temperature: 90°C Detection: DRI Injection: 20 µL

#### Sample:


1. Sucrose 2. Glucose 3. Fructose



#### **Separation of Various Sugars and Sugar** Alcohols on a Coregel-87C Column

#### **Analysis Conditions:**

Column: Coregel-87C (7.8 mm x 300) Eluent: Distilled Water Flow rate: 0.6 mL/min Temperature: 85°C Pressure: 425 psig Detection: RI Range 18x Injection: 20 µL




#### **Analysis of Honey** on a Coregel-87C Column

#### **Analysis Conditions:**

Column: Coregel-87C Eluent: Distilled Water Flow rate: 0.6 mL/min Temperature: 85°C Pressure: 425 psig Detection: RI Range 16x Injection: 20 µL

#### Sample:



3

#### Sugar Separation on CARBOSep CHO-820

#### **Analysis Conditions:** Column: CHO-820 Eluent: Distilled Water Flow rate: 0.5 mL/min Temperature: 90°C Detection: RI 3 Injection: 20mL 1 Sample: 2 1. Melezitose (2.4 mg/mL) 2. Maltose (2.4 mg/mL) 3. Glucose (2.4 mg/mL) 4. Maltitol (3.2 mg/mL) 5. Fucose (2.4 mg/mL) 6. Ribose (2.4 mg/mL) 6 0 5 10 15 20 25

#### **Corn Syrup**

1.DP7


2.DP6

3.DP5

4. DP4


5.DP3


#### **Analysis Conditions:**



#### **Orange Juice**

#### **Analysis Conditions:**






#### **Determination of Sugars in Ale**

#### Analysis Conditions:

Column: CHO-682 Eluent: H<sub>2</sub>O Flow rate: 0.4 mL/min Temperature: 80°C Detection: DRI Injection: 20 µL

#### Sample:

1. Maltose 2. Glucose 3. Ethanol



#### 13


#### Malted Milk Candy

Analysis Conditions:

Column: CHO-411 Eluent:  $H_2O$ Flow rate: 0.5 mL/min Temperature: 75°C Detection: DRI Injection: 20  $\mu$ L of pretreated sample with POLYSorb<sup>TM</sup> ACT-1

#### Sample:

1.DP7 2.DP6 3.DP5 4.DP4 5.DP3 6.Maltose 7.Glucose




#### Mono- and Disaccharides


#### **Analysis Conditions:**

Column: CHO-682 Eluent: H<sub>2</sub>O Flow rate: 0.4 mL/min Temperature: 80°C Detection: Refractometry

#### Sample:

1. Sucrose 2. Maltose 3. Glucose 4. Xylose 5. Galactose 6. Arabinose 7. Mannose






## Saccharides and Sugar Alcohol Separation on CARBOSep CHO-820

#### Analysis Conditions:

Column: CHO-820 Eluent: Distilled Water Flow rate: 0.5 mL/min Temperature: 90°C Detection: RI Injection: 20µL





#### Separation of Carbohydrate Standard

#### **Analysis Conditions:**

Column: CHO-820 Eluent: H<sub>2</sub>O Flow rate: 0.5 mL/min Temperature: 90°C Detection: DRI

#### Sample:





#### CARBOSep CHO-620

#### (6.5 x 300mm) P/N CHO-99-9753

- Calcium form ligand-exchange column
- · Ideal for the separation of monosaccharides and sugar alcohols
- Very reproducible

#### CARBOSep CHO-620 Guard Kit

P/N CHO-99-2353

#### CARBOSep CHO-620 Guard Cartridge – 2/PK

P/N CH0-99-1353

#### CARBOSep CHO-682 Lead

(7.8 x 200mm) P/N CHO-99-8854

#### (7.8 x 300mm) P/N CHO-99-9854

- Lead form ligand-exchange column
- · Ideal for the separation of mono and disaccharides as well as alcohols
- High capacity

#### **CARBOSep CHO-682 Guard Kit**

P/N CHO-99-2354

#### CARBOSep CHO-682 Guard Cartridge – 2/PK

P/N CH0-99-1354

#### CARBOSep CHO-820 Calcium

(7.8 x 200mm) P/N CHO-99-8855

#### (7.8 x 300mm)


- P/N CHO-99-9855
- Calcium form ligand-exchange column
- Designed with balance of resolution and ruggedness


#### **CARBOSep CHO-820 Guard Kit**


P/N CHO-99-2355

## CARBOSep CHO-820 Guard Cartridge – 2/PK

P/N CHO-99-1355







#### CARBOSep CHO-611 OH

#### (6.5 x 150mm) P/N CHO-99-7752

- Sodium form ligand-exchange column
- Designed for use with Sodium Hydroxide eluant
- Compatible with amperometric detection

#### CARBOSep CHO-611 OH Guard Kit

P/N CHO-99-2352

#### CARBOSep CHO-611 OH Guard Cartridge – 2/PK

P/N CH0-99-1352

#### CARBOSep CHO-411

(7.8 x 300mm)

- P/N CHO-99-9850
- Sodium form mixed-mode column
- Separates by both ligand exchange and size exclusion
- Designed for the separation of oligosaccharides up to DP10
- Reproducible separation of corn syrup

#### CARBOSep CHO-611 Guard Kit

P/N CHO-99-2351

#### CARBOSep CHO-611 Guard Cartridge – 2/PK

P/N CH0-99-1351

#### CARBOSep CHO-611

#### (6.5 x 300mm)


- P/N CHO-99-9751
- Sodium form mixed-mode column
- Separates by both ligand exchange and size exclusion
- Designed for the separation of oligosaccharides up to DP5
- Reproducible separation of corn syrup


#### **CARBOSep CHO-611 Guard Kit**


P/N CHO-99-2351

#### CARBOSep CHO-611 Guard Cartridge – 2/PK

P/N CH0-99-1351







ANSGENOMIC BIOCONSUMABL

#### CARBOSep USP L19 CA-FORM

#### (4.0 x 250mm) P/N CHO-99-8453

- Calcium form ligand-exchange column
- Complies with USP L-19 specifications for the separation of sorbitol and mannitol
- Can also separate a wide number of other carbohydrates

#### CARBOSep CHO-820 Guard Kit

P/N CHO-99-2355

#### CARBOSep CHO-820 Guard Cartridge – 2/PK

P/N CHO-99-1355

#### CARBOSep COREGEL-87C

#### (7.8 x 300)

- P/N CHO-99-9860
- Calcium form 9µm ligand exchange resin with 8% cross-linking
- Compatible replacement for the Bio-Rad Aminex HPX 87C
- Designed for the analysis of sugars and sugar alcohols

#### **CARBOSep COREGEL-87C Guard Kit**

P/N CHO-99-2360

#### CARBOSep COREGEL-87C Guard Cartridge – 2/PK

P/N CHO-99-1360

#### **CARBOSep COREGEL-87K**

#### (7.8 x 300)

#### P/N CHO-99-9862


- Potassium form 8µm ligand exchange resin with 8% cross-linking
- Compatible replacement for the Bio-Rad Aminex HPX 87K
- Target application corn syrup and molasses


#### CARBOSep COREGEL-87K Guard Cartridge – 2/PK


P/N CHO-99-1362


#### **Universal Guard Cartridge Holder**

P/N AXC-99-1300









60000-0

#### CARBOSep COREGEL-87N

#### (7.8 x 300mm) P/N CHO-99-9863

- Sodium form 8µm ligand exchange resin with 8% cross-linking
- Compatible replacement for the Bio-Rad Aminex HPX 87N
- Designed for the fast separation of monosaccharides and
- sugar alcohols

#### CARBOSep COREGEL-87N Guard Cartridge – 2/PK

P/N CHO-99-1363

#### **Universal Guard Cartridge Holder**

P/N AXC-99-1300

#### CARBOSep COREGEL-87P

#### (7.8 x 300mm) P/N CHO-99-9864


- Lead form 8µm ligand exchange resin with 8% cross-linking
- Compatible replacement for the Bio-Rad Aminex HPX 87P
- Optimized for the analysis of cellulose hydrolysates


#### CARBOSep COREGEL-87P Guard Cartridge – 2/PK

P/N CHO-99-1364

#### **Universal Guard Cartridge Holder**

P/N AXC-99-1300







#### CARBOSep COREGEL-87MM

#### (7.8 x 300mm) P/N CHO-99-9865

- Mixed calcium/sodium form ligand-exchange column
- Increased efficiency of glucose, fructose, and sugar alcohols
- Easily cleaned with EDTACaNa<sub>2</sub>

#### CARBOSep COREGEL-87MM Guard Cartridge 2/PK

P/N CHO-99-1365


#### **Universal Guard Cartridge Holder**

P/N AXC-99-1300



Eluent: Water Flow rate: 0.6 mL/min Detector: Rl

#### Sample:



## CARBOSep COREGEL-42Ag

#### (7.8 x 300mm) P/N CHO-99-9851

- Silver form ligand-exchange column
- Separate oligosaccharides up to DP11
- Compatible replacement for the Bio-Rad Aminex HPX-42A column

#### CARBOSep COREGEL-42Ag Guard Cartridge 2/PK


P/N CHO-99-1366

#### **Universal Guard Cartridge Holder**

P/N AXC-99-1300

#### Analysis Conditions:

Eluent: H<sub>2</sub>O Flow rate: 0.4 mL/min Detector: RI




# organic ACID Analysis

#### ICSep Columns for Organic Acid Analysis

Ion exclusion is the preferred method for the separation of weakly ionizable species such as organic acids and alcohols. Transgenomic provides a broad range of columns that provide varying efficiencies and selectivities for the separation of weak acids by ion exclusion.

The packings employed with ion exclusion are totally sulfonated polystyrene divinylbenze (PS/DVB) copolymers. By totally sulfonating the polymer, the bead behaves as though it were a negatively charged sphere. This charged sphere is referred to as a Donnan membrane. Species that have a negative charge are repelled from the negatively charged membrane, while uncharged species are allowed to enter the sphere and adsorb onto the beads. The mobile phases employed with ion exclusion are low concentration acids, such as 5mM sulfuric acid.



This equilibrium is regulated by the acidic dissociation constant (pKa) of the organic acid or alcohol. Therefore, species are analyzed by ion exclusion and elute according to their pKa.


#### Features

The key features of the ICSep ion exclusion columns are:

- Polymeric Substrate
- High efficiency
- High resolution
- Separates organic acids, carbohydrates, and alcohols on the same column
- Very Rugged Design which
  provides long life

Since ICSep columns are based on a polymeric substrate consisting of polystyrene/divinylbenzene copolymers they are stable in the pH range of 0 to 14, temperature stable, and very rugged. The ICSep organic acid columns have been shown to last for thousands of runs without cleaning. They show very good lotto-lot and column-to-column reproducibility with retention times varying by less than 1%.

Transgenomic offers ICSep organic acid columns to meet your analytical needs. ICSep columns are available that focus on speed or efficiency and there are ICSep ion exclusion columns that focus on ruggedness and the ability to handle dirty samples. There are even ICSep columns for aromatic organic acids. Transgenomic is sure to have an ion exclusion column to meet your needs.



#### Selectivity Chart for Ion Exclusion Columns


Compound	Coregel 87H @ 85°C	Coregel 64H @ 65°C	ION-300 @ 65°C	ORH-801 @ 45°C	
	(units in minutes)	(units in minutes)	(units in minutes)	(units in minutes)	
Acetic acid	13.8	15.0	14.9	10.4	
Acetoacetic acid	nd	nd	nd	10.2	
Aconitic acid	8.6	9.8	10.7	7.2	
Acrylic acid	15.9	17.7	17.9	13.1	
Adipic acid	12.5	15.1	15.8	11.6	
Butanol	32.9	35.1	25.2	18.4	
Butyric acid	18.4	21.0	20.8	15.2	
Citraconic acid	10.1	11.0	11.5	nd	
Citric acid	7.5	8.0	8.6	5.5	
Ethanol	21.4	21.7	20.6	14.6	
Formic acid	12.9	13.8	13.9	9.6	
Fumaric acid	11.5	13.4	14.7	10.0	
2-Furoic acid	22.1	26.9	29.0	22.0	
Glucoronic acid	nd	nd	nd	5.3	
Glycolic acid	11.4	13.0	12.9	8.5	
Glycoxylic acid	9.2	9.7	10.3	6.5	
Hydroxybutyric acid	12.8	14.0	14.1	9.5	
lso-butyric acid	17.3	19.6	19.5	14.0	
Itaconic acid	11.1	12.8	13.4	9.1	
Keto-butyric acid	nd	nd	11.4	7.4	
Keto-glutaric acid	7.8	8.2	nd	5.6	
Keto-valeric acid	11.7	12.6	13.1	8.6	
Lactic acid	11.9	12.9	11.6	8.7	
Maleic acid	8.2	8.6	9.0	5.9	
Malic acid	8.8	9.6	10.3	6.6	
Malonic acid	9.3	10.0	10.7	6.9	
Methanol	18.7	19.0	18.7	12.9	
Methylglutaric acid	11.8	13.9	14.5	10.0	
Methylsuccinic acid	10.9	12.5	13.0	8.8	
Oxalic acid	6.7	6.6	nd	4.5	
Propanol	25.9	26.7	22.2	16.1	
Propionic acid	15.8	17.4	17.4	12.3	
Pyruvic acid	9.2	9.5	9.9	6.3	
Quinic acid	9.4	10.3	11.4	6.9	
Shikimic acid	10.5	11.8	12.9	8.2	
Succinic acid	10.4	11.7	12.2	8.2	
Tartaric acid	8.0	8.6	9.5	5.9	

Flow rate: 0.6 mL/minute. nd = not determined









#### **Comparison of Organic Acids Retention on Ion-exclusion Columns**

#### Analysis Conditions:

Column: ION-310 (6.5 x 150 mm), 3 ORH-801 (6.5 x 300 mm), ION-300 (7.8 x 300 mm) ION-310 ORH-801 ION-300 Eluent: 0.002 N H<sub>2</sub>SO<sub>4</sub> з Flow rate: 0.5 mL/min Temperature: 35°C Detection: UV at 210 nm 3 Injection: 20 µL Sample: 1. Maleic Acid (2 ppm) 2. Malic Acid (100 ppm) 3. Fumaric Acid (5 ppm) 1 2 2 0 min 5 10 5 10 0 min 15 0 min 5 10 15 20 25


#### **Borate and Bicarbonate**

#### **Analysis Conditions:**

Column: ION-310 Eluent: 0.05 MH<sub>2</sub>SO<sub>4</sub> Flow rate: 0.5 mL/min Temperature: Ambient Detection: Conductivity Injection: 20 µL

#### Sample:

1. Borate (11 ppm) 2. Bicarbonate (60 ppm)




2

#### Wine Analysis by High Resolution **Ion Exclusion**

#### **Analysis Conditions:**

Column: ION-300 Eluent: 0.005 N H<sub>2</sub>SO<sub>4</sub> Flow rate: 0.3 mL/min Temperature: 60°C Detection: DRI

#### Sample:




#### **Analysis of Corn Mash Fermentation Sample**

#### **Analysis Conditions:**

Column: ION-300 Eluent: 0.005 N H<sub>2</sub>SO<sub>4</sub> Flow rate: 0.4 mL/min Temperature: 60°C Detection: UV at 210 Injection: 20 µL filtered

#### Sample:


1. Citric, Isocitric 2. Pyruvic 3. Succinic 4. Fumaric 5. Ethanol



20

24


28



#### **Fast Acid Analysis**

#### Analysis Conditions:

Column: ORH-801 Eluent: 0.01 N H<sub>2</sub>SO<sub>4</sub> Flow rate: 0.5 mL/min Detection: Conductivity




#### **Fluoride in Dental Rinse**

#### **Analysis Conditions:**


Column: ION-310 Eluent: 0.01 N H<sub>2</sub>SO<sub>4</sub> Flow rate: 1.0 mL/min Temperature: 50°C Detection: DRI

#### Sample:

3. Glycerol 5. Ethanol



#### **Separation of Organic Acids**




#### **Determination of Chemical Markers for Thermal Processing of Ground Meat**

1.0xalic

3. Tartaric 4. Malic 5. Lactic

6. Formic 7. Fumaric 8. Propionic 9. Butyric



#### **USP-NF Malic Acid Method**, **Fumaric and Maleic Acids**

#### **Analysis Conditions:**

Column: ORH-801 packing L17 specification Eluent: 0.01 N H<sub>2</sub>SO<sub>4</sub> 2 Flow rate: 0.6 mL/min Temperature: 37°C Detection: UV at 210 nm Injection: 20 µL Sample: USP Malic Acid 3 (100 mg in 100 mL volumetric flask, made up with 0.01 N H <sub>2</sub> SO <sub>4</sub> ) 1. Maleic Acid 2. Malic Acid 3. Fumaric Acid 0 min 10 15 5


#### **Organic Acid Separation** on COREGEL-87H1

#### Analysis Conditions:

Column: COREGEL-87H1 Eluent: 5mM Sulfuric Acid Flow rate: 1.0 mL/min Temperature: 35°C Detection: UV @ 210nm Injection: 20µL

#### Sample:

1. Lactic Acid (2.0 mg/mL) 2. Acetic Acid (2.0 mg/mL) 3. Propionic Acid (2.0 mg/mL) 4. alpha-Butyric Acid (1.0 mg/mL) 5. Glutamic Acid (2.0 mg/mL)



#### **Organic Acid Separation** on COREGEL-87H3

#### Analysis Conditions:


Column: COREGEL-87H3 Eluent: 0.008M Sulfuric Acid Flow rate: 0.6 mL/min Temperature: 35°C Detection: UV @ 210nm Injection: 20µL



**ICSep COREGEL-87H1** 

(7.8 x 100mm) P/N ICE-99-5861

(7.8 x 300mm) P/N ICE-99-9861



#### **ICSep COREGEL 87H Guard Kit**

P/N ICE-99-2361

ICSep COREGEL 87H Guard Cartridge – 2/PK

P/N ICE-99-2371

#### **ICSep ION-300**

#### (7.8 x 300mm) P/N ICE-99-9850


- Select when high resolution is the primary concern
- Separates Organic Acids, Alcohols and Carbohydrates
- all on the same column

#### ICSep GC-801 Guard Kit

P/N ICE-99-2354

#### ICSep GC-801 Guard Cartridge – 2/PK

P/N ICE-99-2364



#### ICSep COREGEL-107H

#### (7.8 x 300mm) P/N ICE-99-9866

- New Higher Cross-linked Column
- Improved Resolution for Organic Acids

ICSep COREGEL-107H Guard Cartridge – 2/PK P/N ICE-99-2366

#### **Universal Guard Cartridge Holder**


P/N AXC-99-1300

#### Organic Acid Separation Comparison on the NEW ICSep COREGEL-107H and Competitive Organic Acid Column

#### Analysis Conditions:

Column: COREGEL-107H and Competitive Organic Acid Column Eluent: 0.008N Sulfuric Acid Flow rate: 0.6 mL/min Temperature: 35°C Detection: UV @ 210nm Injection: 20µL

Sample:
1. Citric Acid
2. Alpha Ketoglutaric Acid
3. Fumaric Acid
4. Acetic Acid



#### **ICSep ORH-801**

#### (6.5 x 300mm) P/N ICE-99-9754


- Provides good balance of high efficiency and ruggedness
- Versatile column for Organic Acids, Alcohols and Carbohydrates

#### **ICSep GC-801 Guard Kit**

P/N ICE-99-2354

#### ICSep GC-801 Guard Cartridge – 2/PK

P/N ICE-99-2364



#### **Sugar and Organic Acid Separation** on ICSep Wine Analysis WA-1

#### **Analysis Conditions:**

Column: Wine Analysis WA-1 Eluent: 0.0025N Sulfuric Acid Flow rate: 0.6 mL/min Temperature: 45°C Detection: RI Injection: 20µL

#### Sample:

1. Citric Acid (0.5 mg/mL) 2. Tartaric Acid (2.0 mg/mL) 3. Glucose (2.0 mg/mL) 4. Malic Acid (1.0 mg/mL) 5. Fructose (2.0 mg/mL) 6. Succinic Acid (0.5 mg/mL) 7. Lactic Acid (2.0 mg/mL) 8. Glycerine (5.0 mg/mL) 9. Acetic Acid (0.5 mg/mL) 10.2,3-Butanediol (0.5 mg/mL) 11. Isomer Impurity 12. Ethanol (10.0 mg/mL)

0

# 10 11

Τ

5

10

| 15

#### **ICSep WA-1 Wine Analysis Column**

(7.8 x 300mm) P/N ICE-99-9810

#### **ICSep WA-1 Wine Guard Kit**

P/N ICE-99-3510

#### **ICSep WA-1 Wine Guard Cartridge 2/PK**

P/N ICE-99-1310

12

#### **ICSep ION-310**

#### (6.5 x 150mm) P/N ICE-99-7752

Designed for fast analysis of organic acids and alcohols


• Ideal for the analysis of borate and bicarbonate

#### ICSep GC-801 Guard Kit

P/N ICE-99-2354

#### ICSep GC-801 Guard Cartridge – 2/PK

P/N ICE-99-2364



#### **ICSep ARH-601**

#### (6.5 x 100mm) P/N ICE-99-5753

• Designed for the separation of Aromatic organic acids

• Uses aqueous mobile phases

#### ICSep GC-601 Guard Kit

P/N ICE-99-2353

#### ICSep GC-601 Guard Cartridge – 2/PK

P/N ICE-99-2363

#### **ICSep COREGEL-64H**


(7.8 x 300mm) P/N ICE-99-9860


#### **ICSep COREGEL 64H Guard Kit**

P/N ICE-99-2360

#### ICSep COREGEL 64H Guard Cartridge – 2/PK

P/N ICE-99-2370





# POLYMERIC REVERSED Phase

#### **RPSep** Columns

Reversed phase is commonly referred to as adsorption chromatography. Reversed phase works by taking advantage of the hydrophobic interactions between molecules and a hydrophobic stationary phase.

In reversed phase, molecules are adsorbed onto a hydrophobic stationary phase. Then, the molecules are desorbed by changing the hydrophobic character of the mobile phase such that the molecules will selectively partition into the mobile phase and elute from the column.


Traditionally, silica-based packings have been the most commonly used sorbants. However, as samples become more challenging, as with biological samples, supports are required that have broader pH ranges, are more rugged, and can be cleaned. Transgenomic provides a family of products all based on polystyrenedivinylbenzene sorbants that utilize our patented alkylation technology.

#### Features

The key features of RPSep polymeric reversed phase columns are:

- pH stable from 0 14
- temperature stable
- very rugged, long lasting materials
- very tight particle size range ( $\pm$  0.5  $\mu m)$  for high efficiency
- very high efficiency for polymeric resins
- both alkylated and non alkylated PS/DVB available
- all resins available in both analytical and bulk for scalability

And, as with all Transgenomic Chromatography products, RPSep columns provide excellent column-to-column and lot-to-lot reproducibility.




#### Aspirin and Salicylic Acid on Poly-RP CO

#### Analysis Conditions:

Column: Poly-RP C0 Eluent: 1% H<sub>3</sub>PO<sub>4</sub> (28%) in 50:50 ACN:H<sub>2</sub>O Flow rate: 0.75 mL/min Temperature: Ambient Detection: UV at 254 nm

#### Sample:

1. Aspirin (2-(acetyloxy)-benzoic acid) 2. Benzoic Acid 3. Salicylic Acid




#### Separation of Sulfonamides on Poly-RP CO

#### Analysis Conditions:

Column: Poly-RP C0 Eluent: 0.01 M KH<sub>2</sub>PO<sub>4</sub> in 25:75 ACN:H<sub>2</sub>O Flow rate: 0.75 mL/min Detection: UV at 254 nm Injection: 10  $\mu$ L

## 1 Sample: 1. Sulfanilic Acid (10 µg/mL) 2. Sulfanilamide (10 µg/mL) 3. Sulfathiazole (20 µg/mL) 4. Sulfamethizole (20 µg/mL) 5. Sulfamethazine (30 µg/mL) 6. Sulfamethazine (30 µg/mL) 7. Sulfisoxazole (30 µg/mL) 8. Sulfamethoxazole (30 µg/mL) 9 Min 5 10 15

#### **Separation of PGRs and Herbicides**




#### **Separation of Triazine Herbicides** on Poly-RP-C0

2. Kinetin

Analysis Conditions:
Column: Poly-RP C0
Eluent: 60:40 ACN:H <sub>2</sub> O
Flow rate: 0.75 mL/min
Temperature: Ambient
Pressure: 107 Bar
Detection: UV at 254 nm

#### Sample: 1. Aminotriazole 2. Simazine 3. Atrazine 4. Propazine 5. Ametryne 6. Prometryne




#### **Carbamates**

#### **Analysis Conditions:**

Column: ACT-1 Eluent: 70:30 ACN:H<sub>2</sub>O Flow rate: 0.5 mL/min Temperature: Ambient Detection: UV at 240 nm Injection: 20 µL

#### Sample:

1. Oxamyl (5 µg/mL) 2. Aldicarb (30 µg/mL) 3. Carbofuran (30 µg/mL) 4. Carbaryl (30 µg/mL) 5. Propham (2.5 µg/mL) 6. Methiocarb (12.5 µg/mL) 7. Ferbam (9 µg/mL) 8. ChlorolPC (9 µg/mL) 9. EPTC (87.5 µg/mL)



# Separation of polar and Non-polar Compounds

#### **Analysis Conditions:** Column: ACT-1 Eluent: 60:40 ACN:H<sub>2</sub>O Flow rate: 0.3 mL/min Temperature: Ambient Detection: UV at 254 nm З 2 2 Sample: 1. Unknown 2. Phenol 3 3. Aniline 4. Acetophenone 5. Nitrobenzene 6.Toluene 0 min 20 30 40 0 min 40 20 30

# **Tertiary Amines on Poly-RP C0**


#### Analysis Conditions:

Column: Poly-RP CO Eluent: 0.1 *M* Ammonia in 80:20 ACN:H<sub>2</sub>O Flow rate: 0.75 mL/min Temperature: Ambient Detection: UV at 210 nm

#### Sample: 0.05 $\mu L/mL$ of 1. Trimethylamine

2. Triethylamine

4. Tripropylamine 5. Tribitylamine




# **Comparison of ACT-1 with PRP-type Column**


## Analysis Conditions:

Column: ACT-1 Eluent: 80:20 Methanol: Water Linear Velocity: 4.2 cm/min Temperature: Ambient Detection: UV at 254 nm

#### Sample:

1. Methylphenone 2. Ethylphenone 3. Propylphenone 4. Butylphenone 5. Pentylphenone





37

# Chromatography

#### Introduction

Ion Chromatography (IC) is the separation of inorganic and organic ionic species by ion exchange chromatography followed by suppressed conductivity detection. The technique was pioneered by Dow Chemical Company in 1974 and has grown in popularity since.

The species analyzed by IC include both anions and cations. The separation of anions is accomplished via anion exchange chromatography. The separations of cations are accomplished via cation exchange chromatography. Transgenomic provides a broad range of columns for the separation of both anions and cations.

The resins used for anion and cation exchange chromatography in IC employ a functionalized, macroporous polystyrene/divinyl benzene copolymer. Resins functionalized with quaternary alkyl or alkynol ammonium groups are used with hydroxide or carbonate-based eluents for anion exchange IC. Resins functionalized with sulfonic acid or carboxylic acid groups are used with acidic eluents for cation exchange IC.

#### Features

The key features of the Transgenomic IC columns are:

- Polymeric substrate
- Solvent compatibility
- High efficiency
- Reproducibility
- pH Stability from 0 to 14

#### **Column Selection**

Transgenomic IC columns have been designed to run on a variety of systems. They are tested to be compatible with Ion Chromatographs from: Metrohm-Peak, Dionex, Hach-Lachat, and Alltech. The selectivities have been optimized to be compatible with many of the common IC columns currently available. This includes columns that meet the requirements of E.P.A. methods 300 parts a and b, and E.P.A. method 300.1.

# **Column Equivalents Guide**


TRANSGENOMIC COLUMN	COMPETITIVE COLUMNS	APPLICATION
ICSep AN300	Dionex AS4A	F <sup>-</sup> , Cl <sup>-</sup> , NO <sub>2</sub> <sup>-</sup> , Br <sup>-</sup> , NO <sub>3</sub> <sup>-</sup> , HPO <sub>4</sub> <sup>2<sup>-</sup>, SO<sub>4</sub><sup>2<sup>-</sup>, By E.P.A. Method 300.0(a)</sup></sup>
ICSep AN1	Dionex AS9-HC	$F^-$ , $CI^-$ , $NO_2^-$ , $Br^-$ , $NO_3^-$ , $HPO_4^{2^-}$ , $SO_4^{2^-}$ , Low molecular weight, Organic acids in medium to high ionic strength matrices
		Cr(III), Cr(VI) as CrO <sub>3</sub> -, CrO <sub>4</sub> <sup>2-</sup>
ICSep ANSC	Dionex AS4A-SC	Polyvalent Phosphates, Arsentate, Sulfite Selenate, Arsenite, Selenite, $F^-$ , $CI^-$ , $NO_2^-$ , $Br^-$ , $NO_3^-$ , $HPO_4^{2^-}$ , $SO_4^{2^-}$ , Low molecular weight, Organic acids
ICSep AN1SC	Dionex AS9-HC	$F^-$ , $CI^-$ , $NO_2^-$ , $Br^-$ , $NO_3^-$ , $HPO_4^{2^-}$ , $SO_4^{2^-}$ , Low molecular weight, Organic acids in medium to high ionic strength matrices
ICSep AN2	Dionex AS14	Arsenate, Sulfite, Selenate, Arsenite, Selenite $F^-$ , $CI^-$ , $NO_2^-$ , $Br^-$ , $NO_3^-$ , $HPO_4^{2^-}$ , $SO_4^{2^-}$ , Low molecular weight Organic acids
ICSep AN300B	Dionex AS9	F <sup>-</sup> , Cl <sup>-</sup> , NO <sub>2</sub> <sup>-</sup> , Br, NO <sub>3</sub> <sup>-</sup> , HPO <sub>4</sub> <sup>2-</sup> , SO <sub>4</sub> <sup>2-</sup> , ClO <sub>2</sub> <sup>-</sup> , ClO <sub>3</sub> <sup>-</sup> , BrO <sub>3</sub> <sup>-</sup>
ICSep CN2	Dionex CS15	Li+, Na+, K+, Rb+, Cs+, Mg <sup>2+</sup> , Ca <sup>2+</sup> , NH <sup>4+</sup> , Cu <sup>2+</sup> , Ni <sup>2+</sup> , Zn <sup>2+</sup> , Co <sup>2+</sup> , Cd <sup>2+</sup> , Pb <sup>2+</sup> , Mn <sup>2+</sup> , Fe <sup>2+</sup> , Fe <sup>3+</sup>

# Anions by E.P.A. Method 300.0(a)

#### Conditions

Column: ICSep AN300 Eluent: 1.7mM Sodium Carbonate, 1.8mM Sodium Bicarbonate Flow rate: 2.0 mL/min Detection: suppressed conductivity

#### Sample:






# Anions by E.P.A. Method 300.1

#### Conditions

Column: ICSep AN300B Eluent: 3.5mM Sodium Carbonate Flow rate: 1.0 mL/min Detection: conductivity

#### Sample:




# **Anion Separation using ICSep ANSC**

#### Conditions

Column: ICSep ANSC Eluent: 1.8mM Sodium Carbonate, 1.7mM Sodium Bicarbonate Flow rate: 1.2 mL/min Detection: suppressed conductivity

#### Sample:




#### **Determination of Perchlorate** using ICSep ANSC

#### Conditions

Column: ICSep ANSC with guard Eluent: 30mM Sodium Hydroxide, 10mM Cyanophenol Flow rate: 1.2 mL/min Detection: suppressed conductivity





# **Cations using ICSep CN2**


#### Conditions

Column: ICSep CN2 Eluent: 0.1mM Ce (III) Flow rate: 1.0 mL/min Detection: UV @ 254nm

#### Sample:

1.3ppm sodium 2.3ppm ammonium 3.5ppm potassium 4.30ppm rubidium 5.30ppm cesium

7.10ppm calcium



# **Ordering Information**

DESCRIPTION	PART NUMBER
ICSep AN2, 4.6mm x 250mm	ANX-99-8515
ICSep AN2 Guard Column, 4.6mm x 50mm	ANX-99-3515
ICSep AN2 Guard Cartridges, 3/pk, 3.0mm x 10mm	ANX-99-0015
ICSep AN1, 4.6mm x 250mm	ANX-99-8511
ICSep AN1 Guard Column, 4.6mm x 50mm	ANX-99-3510
ICSep AN1 Guard Cartridges, 3/pk, 3.0mm x 10mm	ANX-99-0010
ICSep AN1-SC, 4.6mm x 250mm	ANX-99-8514
ICSep AN1-SC Guard Column, 4.6mm x 50mm	ANX-99-3514
ICSep AN1-SC Guard Cartridges, 3/pk, 3.0mm x 10mm	ANX-99-0014
ICSep AN300, 5.5mm x 150mm	ANX-99-7613
ICSep AN1 Guard Column, 4.6mm x 50mm	ANX-99-3510
ICSep AN1 Guard Cartridges, 3/pk, 3.0mm x 10mm	ANX-99-0010
ICSep AN300B, 4.6mm x 250mm	ANX-99-8516
ICSep AN300B Guard Column, 4.6mm x 50mm	ANX-99-3516
ICSep AN300B Guard Cartridges, 3/pk, 3.0mm x 10mm	ANX-99-0016
ICSep ANSC, 4.6mm x 250mm	ANX-99-8512
ICSep ANSC Guard Column, 4.6mm x 50mm	ANX-99-3512
ICSep ANSC Guard Cartridges, 3/pk, 3.0mm x 10mm	ANX-99-0012
ICSep ION-120, 4.6mm x 120mm	ANX-99-6550
ICSep ION-120 Guard Kit, 4.0mm x 24mm	ANX-99-2350
ICSep ION-120 Guard Cartridges, 3/pk, 4.0mm x 24mm	ANX-99-0090
ICSep CN2, 3.2mm x 100mm	CTX-99-5250
ICSep CN2 Guard Cartridges, 2/pk, 3.0mm x 10mm	CTX-99-1350
ICSep CN2 FA, 4.6mm x 50mm	CTX-99-3550
ICSep CN2 Guard Cartridges, 2/pk, 3.0mm x 10mm	CTX-99-1350

# GUARD-DISC® PROTECTION System

#### **Guard-Disc System**

The Transgenomic Guard-Disc System is a patented column protection system that is designed to provide the protection capabilities of a guard column without adding any extra volume that might interfere with chromatographic separation.

The Guard-Disc System is comprised of a disc, which is available in a variety of functionalities, and a disc holder that couples directly to the column.

The disc is a PEEK ring that contains a functionalized chromatographic membrane. This chromatographic membrane is available in a variety of stationary phases for both HPLC and Ion Chromatography applications.

#### **Phases**


The stationary phases that Guard-Discs Systems are available in include:

- C18
- C8
- Styrene/DVB
- Anion Exchange
- Cation Exchange

It is these functional groups that bind the contaminants that would otherwise be trapped on your analytical column.

#### **Double Protection**

Transgenomic Guard-Disc Systems are porous as well. Not only do they bind species that may contaminate your analytical column, they also filter out particulates that would otherwise be trapped on your analytical column. The Transgenomic Guard-Disc System provides double protection for your chromatographic column.



Membrane Functionality	Application	Porosity (µm)	Solvent Compatibility	pH Range
C18-A	Reversed Phase	0.2	All	2-8
C18-B	Reversed Phase	0.8	Acetonitrile Methanol	2-8
C8	Reversed Phase	0.2	All	2-8
S/DVB	Reversed Phase	0.2	All	1-13
ANEX	Anion Exchange	0.2	All	1-13
CATEX	Anion Exchange	0.2	All	1-13

# **Guard-Disc System Characteristics**

# TRANSGENOMIC GUARD Discs®

Ion Exchangers

# ANEX Guard-Disc – (10/pk)

P/N GRD-99-0704 CATEX Guard-Disc – (10/pk)

P/N GRD-99-0705

## Adsorbants

C18A Guard-Disc (10/pk)

#### P/N GRD-99-0701

C18B Guard-Disc (10/pk) P/N GRD-99-0731

# C8 Guard-Disc (10/pk)

P/N GRD-99-0702

S/DVB Guard-Disc (10/pk)

P/N GRD-99-0706

# TRANSGENOMIC GUARD Disc<sup>®</sup> Holders

# **Guard-Disc Direct Holder 1**

(Parker Type) P/N AXC-99-0002

# **Guard-Disc Direct Holder 2**

(Waters Type) P/N AXC-99-0003

Guard-Disc Universal Holder 1N

(Universal) P/N AXC-99-0004

# solid phase Extraction

# **Transgenomic POLYSorb<sup>™</sup> Products** for Solid Phase Extraction

Solid Phase Extraction (SPE) is a sample preparation technique that is employed to clean up or concentrate samples prior to analysis. SPE can be used to clean-up samples by removing interferences that would otherwise compromise analysis. It can be used to concentrate by allowing a large volume of sample to be reduced into a small elution volume. Compared to other sample preparation techniques, such as liquid-liquid extraction, SPE provides cleaner extracts with high recoveries. SPE is also faster and uses less solvent which saves money.

#### Modes

SPE tubes can be used in two modes:

- In the flow-through mode the sample can be passed through the tube. While passing through the tube, the contaminants present are retained while the analyte of interest is allowed to pass through. The steps for this mode are 1) Load the sample into the tube
   Wash to elute the analyte of interest.
- In the selective elution mode the sample is passed through the tube. But in this mode, the analyte of interest is retained while contaminants pass through. After the sample is loaded onto the column, the analyte of interest is selectively eluted by choosing elution conditions that will elute the analyte from the column while retaining interfering components. The steps used with this mode are 1) Load the sample onto the column 2) Wash through weakly retained or unretained contaminants 3) Elute the analyte of interest.

The most common SPE packing are polar adsorbants. These adsorbants are used to remove organic interferences from samples. Also, commonly used are ion exchangers to remove charged species as interferences. Transgenomic offers products for both adsorption and ion exchange.

#### Key Features of Transgenomic SPE products

As with all of Transgenomic's chromatography products, the SPE products are all based on polymeric resins. Polymer-based resins are used because of the broad pH range available and the chemical and physical stability of the materials. These cartridges are ideally suited for cleaning up samples in tough matrices.

Transgenomic POLYSorb cartridges provide very high loading capacities to accommodate for concentrated samples. POLYSorb cartridges also provide excellent selectivity even for trace level analysis.

#### POLYSorb Cartridges in the format you need

Transgenomic POLYSorb cartridges are provided in three stationary phase formats:

- Unmodified Poly-[styrene/divinylbenze] (PS/DVB)
- Alkylated (C18) PS/DVB
- Sulfonated PS/DVB

Transgenomic offers each of these cartridges in either 100mg or 400mg tubes, or we can custom pack in sizes to meet your specific needs.

POLYSorb tubes are compatible with off-the-shelf SPE vacuum manifolds, automated workstations or other commonly used accessories.

# **Extraction of Organic Acids from Burgundy Wine with ACT-1**

#### **Sample Preparation:**

Dilute wine 1:10 with distilled water

#### **Conditioning Step:**

Wet tube with 1 mL of methonal followed by 1 mL of 10:90 methonal:water

Sample Addition:

Load 500 µL of dilute wine

Wash Step: 1.0 mL of water

**Elution Step:** 


1.0 mL of 0.05 N H<sub>2</sub>SO<sub>4</sub>

#### **Analysis Conditions:**

Column: ION-300 Eluent: 0.01 N H<sub>2</sub>SO<sub>4</sub> Flow rate: 0.5 mL/min Temperature: 60°C Detection: UV at 214 nm Injection: 20 µL

#### Sample:

1. Citric Acid 2. Tartaric Acid 3. Glucose 4. Malic Acid 5. Fructose 6. Glycerol 7. Succinic Acid 8. Acetic Acid



Before Sample Prep



45

# POLYSorb ACT-1, C18, 100mg

(100/box) P/N SPE-99-0100

# POLYSorb ACT-1, C18, 400mg

#### (50/box) P/N SPE-99-0101

- Patented, Octadecyl-Alklyated PS/DVB
- Ideal for removal of polar compounds
- Stable over pH 0-14, very rugged

# POLYSorb, MP-3, Highly Sulfonate, 100mg

(100/box) P/N SPE-99-0104

# POLYSorb, MP-3, Highly Sulfonated, 400mg

# (50/box)

P/N SPE-99-0105

- pH stable cation exchange resin
- · Ideal for removing amines
- Remove cations from ICP analysis

# POLYSorb, MP-DVB, PS/DVB 100mg

(100/box)

P/N SPE-99-0108

# POLYSorb, MP-DVB, PS/DVB 400mg

#### (50/box)

# P/N SPE-99-0109

- Non-functionalized styrene-divinylbenzene
- Ideal for removing polar compounds
- pH stable from 0-14
- Also available in bulk

# BULK POLYMERIC Resin

Transgenomic has scale-up in mind every time we develop a new resin. The resin in any column discussed in this catalogue is also available in bulk. This allows you to pack your own analytical columns, then quickly and easily scale your analytical application to semi-prep and preparative scales without redevelopment.

If we do not have the resin or particle size that you need, simply call. We have over 20 years experience in the development of polymer materials for analytical and preparative chromatography applications; allow us to put our expertise to work for you.

# BUFFERS and SOLVENTS FOR HPLC

# **Buffers and Solvents for Reversed Phase Chromatography**

Part Number	Description	Size
56011	Acetonitrile, HPLC Grade	1 liter
700002	Water, HPLC Grade	4 liter
553303	Triethlammonium acetate solution, 2M	200 mL
SP5890	Triethlammonium acetate solution, 2M	6 x 200 mL

# **Amino Acid Analysis Buffers**

Part Number	Description	Size
AAA-99-4086	Sodium Diluent Na200	4 liter
AAA-99-4081	Sodium Eluent Na315	4 liter
AAA-99-4096	Sodium Eluent Na740	4 liter
AAA-99-4085	Sodium Regenerant RG011	4 liter

Custom Amino Acid Buffers are available for your analysis, please contact Transgenomic for further information

47

# HPLC COLUMN Hardware

# **Column Coupler**

The patented Column Coupler was developed for the demanding constraints of high efficiency HPLC columns. The Column Coupler permits the quick and easy connection of two analytical HPLC columns in series, or direct connection between a Valco injection valve and an analytical column. Seals are rated to 5,000psi

The unit is a precision-machined, double-ended PEEK connector with 10-32 threads and a non-wetted Delrin® knurled body. The inert composition and the large knurled handle allow easy, fingertight connections and leakproof seal to 5,000psi. The 0.010" through-hole minimizes extra column volume effects and is compatible with the demanding constraints imposed with use of 3µm packing and microbore HPLC. These couplers are not capable of universal applications since the tip sizes are fixed



# **Guard Cartridge Holder**

The Universal Guard Cartridge Holder was designed for use with Transgenomic guard cartridges.

# Ordering Information:

Part Number	Description
282013	Column Coupler, PEEK
AXC-99-1300	Universal Guard Cartridge Holder, 4.0mm x 24mm

The unit is a stainless steel body with dimensions of 4.0mm x 24mm



Amino Acid Columns	PAGE
Transgenomic Na +	6
Transgenomic Li +	6
Transgneomic Na + Column for System Gold	6
AMINOSep AA-911	7
AMINOSep AA-511	7
AMINOSep AA-511High Speed	7

Carbohydrate Analysis Columns	PAGE
CARBOSep CHO-620	16
CARBOSep CHO-682	16
CARBOSep CHO-820	16
CARBOSep CHO-611 OH	17
CARBOSep CHO-611	17
CARBOSep CHO-411	17
CARBOSep USP-L19	18
CARBOSep COREGEL-87C	18
CARBOSep COREGEL-87K	18
CARBOSep COREGEL-87N	19
CARBOSep COREGEL-87P	19
CARBOSep COREGEL-87MM	20
CARBOSep COREGEL-42Ag	20

Organic Acid Analysis Columns	PAGE
ICSep COREGEL-87H	27
ICSep ION-300	28
ICSep COREGEL-107H	29
ICSep ORH-801	30
ICSep WA-1 Wine Analysis	30
ICSep ION-310	31
ICSep ARH-601	31
ICSep COREGEL-64H	31

Polymeric Reversed Phase Columns	PAGE
RPSep POLYRP CO	37
RPSep ACT-1	37
RPSep PRX-1	37

Ion Analysis Columns	PAGE
ICSep AN1	41
ICSep AN1S	41
ICSep AN2	41
ICSep ANSC	41
ICSep AN300	41
ICSep AN300B	41
ICSep ION-120	41
ICSep CN2	41





# TRANSGENOMIC<sup>®</sup> BIOCONSUMABLES<sup>®</sup>

#### To order catalog items:

TRANSGENOMIC CORPORATE HEADQUARTERS 12325 Emmet Street Omaha, NE, USA 68164

#### Tel: 1-888-233-9283 Fax: 1-402-452-5401 E-mail: orders@transgenomic.com

#### For technical support:

TRANSGENOMIC – SAN JOSE OFFICE 2032 Concourse Drive San Jose, CA 95131

Tel: 1-800-605-0267 Fax: 1-408-432-3231

www.transgenomic.com

