

ACE® Analytical Columns 產品規格表

Bonded Phase	USP Listing	Particle Size (µm)	Pore Size (Å)	pH Range	Temp. Max °C (low pH)	Temp. Max °C (high pH)	Carbon Load (%)	Surface Area (m ² /g)
C18	L1	3 5 10	100	2-8	60	40	15.5	300
C18-AR	L1	3 5 10	100	2-8	60	40	15.5	300
C18-PFP	L1	3 5 10	100	2-8	60	40	14.3	300
C18-HL	L1	3 5 10	90	2-8	60	40	20	400
C8	L7	3 5 10	100	2-8	60	40	9	300
C4	L26	3 5 10	100	2-8	60	40	5.5	300
Phenyl	L11	3 5 10	100	2-8	60	40	9.5	300
CN	L10	3 5 10	100	2-7	60	40	5.5	300
AQ	L1	3 5 10	100	2-8	60	40	14	300
SIL	L3	3 5 10	100	2-7	40	40	N/A	300