

Excel UHPLC Columns 產品規格表

Bonded Phase	USP Listing	Particle Size (µm)	Pore Size (Å)	pH Range	Temp. Max °C (low pH)	Temp. Max °C (high pH)	Carbon Load (%)	Surface Area (m <sup>2</sup> /g)	Endcapped
C18	L1	1.7 2 3 5	100	2-8	60	40	15.5	300	YES
SuperC18	L1	1.7 2 3 5	90	1.5-11.5	60	40	14.8	400	Encapsulated
C18-AR	L1	1.7 2 3 5	100	2-8	60	40	15.5	300	YES
C18-PFP	L1	1.7 2 3 5	100	2-8	60	40	14.3	300	YES
C18-Amide	L1	1.7 2 3 5	100	2-8	60	40	16.4	300	YES
CN-ES	L1	1.7 2 3 5	100	2-8	60	40	12.6	300	YES
C8	L7	2 3 5	100	2-8	60	40	9	300	YES
C4	L26	2 3 5	100	2-8	60	40	5.5	300	YES
Phenyl	L11	2 3 5	100	2-8	60	40	9.5	300	YES
CN	L10	2 3 5	100	2-7	60	40	5.5	300	YES
AQ	L1	2 3 5	100	2-8	60	40	14	300	YES
NH <sup>2</sup>	L8	1.7 3 5	100	2-7	60	40	3.5	300	YES
SIL	L3	2 3 5	100	2-7	40	40	N/A	300	YES