

Adsorbosphere 系列 Dr. Maisch 產品規格表

Phase	Particle Size (µm)	Pore Size (Å)	Surface Area (m ² /g)	Carbon Loading	Endcapping	USP
C18	3 5 10	80	200	12%	YES	L1
C18 HS	3 5	60	350	20%	YES	L1
C18 UHS	5 10	60	500	30%	YES	L1
C8	5	80	200	8%	YES	L7
Phenyl	5	80	200	5%	YES	L11
Cyano	5	80	200	-	YES	L10
Cyano-AQ	5	120	170	-	NO	L10
Amino (NH ₂)	3 5	80	200	-	NO	L8
Silica	5	80	200	-	-	L3
SAX	5	80	200	-	-	-
XL C18	3 5	90	200	11%	YES	L1
XL C18-B	5	90	200	12%	YES	L1
XL C8	3 5	90	200	6%	YES	L7
XL C1 (TMS)	5	90	200	-	YES	L13
XL Silica	5	90	200	-	-	L3
XL SAX	5 10	90	200	-	YES	-
XL SCX	5 10	90	200	-	YES	-