

Agilent GC 管柱 產品規格表

PHASE	填充材	SIMILAR PHASE
DB-1	100% Dimethylpolysiloxane	HP-1, HP-101, Ultra-1, SPB-1, OV-1, 007-1, SP-2100, Rtx-1, XP-1
DB-5	(5%-Phenyl)-methylpolysiloxane	HP-5, Ultra-2, SPB-5 CP-Sil 8CB, RSL-200, Rtx-5 BP-5, OV-5, XP-5
DB-5ms	Equivalent to a (5%-Phenyl)-methylpolysiloxane	Ultra-2, Rtx-5MS, HP-5MS, PTE-5
DB-17	(50%-Phenyl)-methylpolysiloxane	HP-17, RSL-300, Rtx-50, PE-17, 007-17 (MPS-50), SP-2250, SPB-50
DB-17ht	(50%-Phenyl)-methylpolysiloxane	Rtx-65TG
DB-23	(50%-Cyanopropyl)-methylpolysiloxane	SP-2330, Rtx-2330, 007-23
DB-35	(35%-Phenyl)-methylpolysiloxane	Rtx-35, SPB-35, AT-35, Sup-Herb
DB-35ms	Equivalent to a (35%-Phenyl)-methylpolysiloxane	Rtx-35, SPB-35, AT-35, Sup-Herb
DB-200	(35%-Trifluoropropyl)-methylpolysiloxane	RSL-400, SP-2401, Rtx-200
DB-225	(50%-Cyanopropylphenyl)-methylpolysiloxane	HP-225, SP-2330, CP-Sil 43CB, RSL-500, Rtx-225, OV-225, XP-225
DB-502.2		Rtx-502.2, VOCOL, HP-VOC
DB-5.625	Equivalent to a (5%-Phenyl)-methylpolysiloxane	Ultra-2, XTI-5, Rtx-5, PTE-5, HP-5Trace
DB-608		SPB-608, NON-PAKD Pesticide, PE-608, 007-608, HP-608
DB-624	mid-polar 6% cyanopropyl/phenyl 94% polydimethylsiloxane	VOCOL, NON-PAKD AT-624, HP-624, Rtx-Volatiles, XP-624
DB-Dioxin		SP-2331, 007-23, Rtx-2332, CP Sil 88

Agilent GC 管柱 產品規格表

PHASE	填充材	SIMILAR PHASE
DB-1301	(6%-Cyanopropylphenyl)-methylpolysiloxane	Rtx-1301, CP-624, PE-1301, HP-1301
DB-1701	(14%-Cyanopropylphenyl)-methylpolysiloxane	SPB-1701, CP-Sil 19CB, Rtx-1701, OV-1701, PE-1701, XP-1701
DB-EVDX		Equivalent to (5%-Phenyl)-methylpolysiloxane
DB-FFAP	Nitroterephthalic Acid Modified Polyethylene Glycol	Stabilwax DA, XP-FFAP
DB-Petro	100% Dimethylpolysiloxane	Petrocol DH, HP-PONA, SPB-1, 007-1, Rtx-1, MXT-1
DB-TPH		
DB-2887	100% Dimethylpolysiloxane	HP-1, Petrocol EX2887, MXT-1, MXT-2887, HT Sim Dist SB
DB-VRX		VOCOL, NON-PAKD AT-624, HP-624, Rtx-Volatiles, XP-624
DB-XLB		Rtx-XLB, MDN-12, ZB-XLB, ZB-XLB HT
DB-WAX	Polyethylene Glycol (PEG)	HP-20M, SUPELCOWAX 10, Stabilwax, HP-Innowax, XP-WAX
DB-WAXetr	Polyethylene Glycol (PEG)	HP-20M, SUPELCOWAX 10, PE-CW, CP-WAX 52CB, HP-Innowax
Cyclodex-B	10.5% β -cyclodextrin in DB-1701	LIPODEX C, Rt- β DEXm, β -DEX 110, β -DEX 120
CAM	Base Deactivated Polyethylene Glycol	Stabilwax-DB, Carbowax Amine, HP-20M
GS-Alumina (PLOT)	Alumina	Al ₂ O ₃ /KCl, Al ₂ O ₃ /Na ₂ SO ₄ , Rt-Alumina PLOT, Alumina PLOT