

Synthetic adsorbents and reversed phase media 產品規格表

Product Name	Base Material	Mean particle size (μm)	Pore Size	pH range	Typical Application
CHP85 / P120	St / DVB	120	Middle	Whole range	Small molecules extraction
CHP87 / P120	St / DVB	120	Small	Whole range	Small molecules extraction
CHPOD / P30	St / DVB	30	Large	2-12	Small molecules extraction
CSP800	St / DVB	120	Middle	Whole range	Adsorption of non-ionic substances
SFP08 / P25	St / DVB	25	Middle	Whole range	Small molecules extraction