

## GLD Pro Gas Leak Detector

Aids in quickly locating and identifying gas leaks

The Thermo Scientific™ GLD Pro™ Gas Leak Detector is specifically designed for use with gas chromatography instruments. Detection of leaks allows the user to reduce detector noise, provide a stable baseline, reduce carrier gas by minimizing waste, and maximize the lifetime of the analytical column by minimizing the presence of oxygen and other impurities in the carrier gas.

- Suitable for detection of a wide range of laboratory gases
- Push button on/off switch
- Push button zero function
- Automatic shutoff (5 minutes)
- LED light indicator for intensity of leak
- Rechargeable battery (up to 6 hrs. operation)
- Durable storage case
- Probe holder
- One year warranty



### GLD Pro Gas Leak Detector

Description	Cat. No.	Quantity
GLD Pro Leak Detector	<b>66002-001</b>	1 Each
Small Probe Adaptor	<b>66002-003</b>	1 Each
Soft-Sided Carry Case (Leak Detector not included)	<b>66002-002</b>	1 Each

