

GC Packed columns 充填分離管柱

固相支撐體 (Solid Support) 之選擇 :

SUPPORT	MESH
Chromosorb W-HP	60/80, 80/100, 100/120, 120/140
Chromosorb G-HP	60/80, 80/100, 100/120
Chromosorb 750	60/80, 80/100, 100/120
Chromosorb 101, 102, 103, 104, 105, 106, 107, 108	60/80, 80/100, 100/120
Chromosorb P (NAW), (AW)	60/80, 80/100, 100/120
Chromosorb G (NAW), (AW)	60/80, 80/100, 100/120
Chromosorb W (NAW), (AW)	60/80, 80/100, 100/120
Chromosorb P (HMDS), (AW-DMCS)	60/80, 80/100, 100/120
Chromosorb G (AW-DMCS)	60/80, 80/100, 100/120
Chromosorb W (HMDS), (AW-DMCS)	60/80, 80/100, 100/120
Chromosorb T	30/60, 40/60
Gas Chrom Q	60/80, 80/100, 100/120, 120/140
Gas Chrom P, Z, A, RA, RZ	60/80, 80/100, 100/120
Porasils, Porapak	50/80, 80/100

各廠牌固相支撐體 (Solid Support) 材質對照表

Equivalent Supports				
MANUFACTURER	NON-ACID WASHED	ACID WASHED	ACID WASHED DMCS TREATED	ACID WASHED DMCS TREATED HIGH PERFORMANCE
Analabs		Anakrom A	Anakrom AS	Anakrom Q
Applied Science	Gas Chrom S, Gas Chrom R	Gas Chrom A, Gas Chrom RA	Gas Chrom Z, Gas Chrom RZ	Gas Chrom Q
Johns Manville	Chromosorb W-NAM	Chromosorb W-AW	Chromosorb W-AW DMCS	Chromosorb W-HP
Johns Manville	Chromosorb G-NAW	Chromosorb G-AW	Chromosorb G-AW DMCS	Chromosorb G-HP
Supelco				Supelcoport