

GC Packed columns 充填分離管柱

其他特殊充填材料

MOLECULAR SIEVES, ALUMINA, AND SILICA GEL		
CAT. NO.	DESCRIPTION	QTY
MOLECULAR SIEVE 4A		
977102	40/60 mesh	100g
977103	60/80 mesh	100g
977104	80/100 mesh	100g
977105	100/120 mesh	100g
MOLECULAR SIEVE 5A		
977106	40/60 mesh	100g
977107	60/80 mesh	100g
977108	80/100 mesh	100g
977109	100/120 mesh	100g
MOLECULAR SIEVE 13X		
977110	40/60 mesh	100g
977111	60/80 mesh	100g
977112	80/100 mesh	100g
977113	100/120 mesh	100g
ACTIVATED ALUMINA (Alcoa Type F-1)		
977118	40/60 mesh	100g
977119	60/80 mesh	100g
977120	80/100 mesh	100g
977121	100/120 mesh	100g
SILICA GEL (Davison Grade 12)		
977122	40/60 mesh	100g
977123	60/80 mesh	100g
977124	80/100 mesh	100g
977125	100/120 mesh	100g

POROUS POLYMERS PORAPAK®			
Porapak is supplied in packs sufficient to fill 10 column 6' x 8" O.D.			
PORAPAK	CAN. NO		QTY
	80/100 MESH	100/120 MESH	
P	977258	977266	20gms
P-S	977259	977267	20gms
Q	977260	977268	26gms
Q-S	977261	977269	26gms
R	977262	977270	24gms
S	977263	977271	26gms
N	977264	977272	29gms
T	977265	977273	31gms

POROUS POLYMERS HAYESEP®			
TYPE	MAX TEMP	CAT. NO	MESH
P	250°C	900003	60/80
		900004	80/100
		900005	100/120
Q	275°C	900000	60/80
		900001	80/100
		900002	100/120
R	250°C	900006	60/80
		900007	80/100
		900008	100/120
S	250°C	900009	60/80
		900010	80/100
		900011	100/120
N	165°C	900015	60/80
		900016	80/100
		900017	100/120
T	165°C	900012	60/80
		900013	80/100
		900014	100/120
A	165°C	900018	60/80
		900019	80/100
		900020	100/120
B	190°C	900022	60/80
		900023	80/100
		900024	100/120
C	250°C	900025	60/80
		900026	80/100
		900027	100/120
D	250°C	900028	60/80
		900029	80/100
		900030	100/120
DB	—	900037	60/80
		900038	80/100
		900039	100/120

TENAX®

Based on 2, 6-diphenyl-p-phenylene oxide, TENAX is suitable for the separation of high boiling compounds such as alcohols, diols, polyethylene glycol compounds, phenols, mono and diamines, ethanalamines, amides, aldehydes, and ketones.
Maximum operating temperature: 350°C

CAT. NO	DESCRIPTION	QTY
977130	TENAX 35/60 mesh	15g
977131	TENAX 60/80 mesh	15g
977132	TENAX 80/100 mesh	15g