

GC Packed columns 充填分離管柱

GC 充填分離管柱快速查閱指南

以下所列充填分離管柱為市面常見材質之充填分離管，
所列 CAT NO 為 6' x 1/8"不鏽鋼材質規格，
供您快速查閱！

STOCK STAINLESS STEEL GC COLUMNS		This is a general configuration. The 1/8" tubing can be carefully bent to fit your GC.
DESCRIPTION	CAT. NO.	
10% Apizon L on Chromosorb W-HP 80/100	975600	
10% Carbowax 400 on Chromosorb W-HP 80/100	975602	
10% Carbowax 1500 on Chromosorb W-HP 80/100	975604	
10% Carbowax 20 M on Chromosorb W-HP 80/100	975606	
10% Carbowax 20 M on Chromosorb W-AW 80/100	975607	
10% Carbowax 20 M-TPA on Chromosorb W-HP 80/100	975608	
Chromosorb 101, 80/100	975610	
Chromosorb 102, 80/100	975612	
Chromosorb 103, 80/100	975614	
Chromosorb 104, 80/100	975616	
Chromosorb 105, 80/100	975618	
Chromosorb 106, 80/100	975620	
Chromosorb 107, 80/100	975622	
Chromosorb 108, 80/100	975624	
10% DEGS on Chromosorb W-AM 80/100	975626	
15% DEGS on Chromosorb W-AM 80/100	975628	
5% FFAP on Chromosorb W-HP 80/100	975630	
5% FFAP on Chromosorb W-HP 80/100	975632	
3% OV-1 on Gas Chrom Q 80/100	975633	
3% OV-1 on Chromosorb W-HP 80/100	975634	
3% OV-1 on Chromosorb W-HP 100/120	975636	
10% OV-1 on Chromosorb W-HP 80/100	975638	
3% OV-17 on Chromosorb W-HP 80/100	975640	
3% OV-17 on Chromosorb W-HP 100/120	975642	
3% OV-17 on Gas Chrom Q 80/100	975644	
3% OV-17 on Gas Chrom Q 100/120	975645	
10% OV-17 on Chromosorb W-HP 80/100	975646	
10% OV-17 on Chromosorb W-HP 100/120	975648	
3% OV-101 on Gas Chrom Q 80/100	975649	
3% OV-101 on Chromosorb W-HP 80/100	975650	
3% OV-101 on Chromosorb W-HP 100/120	975652	
10% OV-101 on Chromosorb W-HP 80/100	975654	
3% OV-210 on Chromosorb W-HP 80/100	975656	
3% OV-225 on Chromosorb W-HP 100/120	975658	
3% OV-225 on Chromosorb W-HP 80/100	975659	
10% OV-225 on Chromosorb W-HP 80/100	975684	
3% OV-275 on Chromosorb W-AW 80/100	975660	
3% SE-30 GC Grade on Chromosorb W-HP 80/100	975661	
Porapak Q, 80/100	975662	
Porapak P, 80/100	975664	
Porapak R, 80/100	975666	
Porapak S, 80/100	975668	
Porapak T, 80/100	975670	
Porapak T, 100/120	975672	
Porapak N, 80/100	975674	
Porapak QS, 80/100	975676	
Porapak PS, 80/100	975678	
1.5% OV-17 + 1.95% QF-1 on Chromosorb W-HP 80/100	975680	
4% SE-30 + 6% OV-210 on Chromosorb W-HP 80/100	975682	

GC Packed columns 充填分離管柱

GC 充填分離管柱快速查閱指南

以下所列充填分離管柱為市面常見材質之充填分離管，
所列 CAT NO 為 6' x 1/8" 不鏽鋼材質規格，
供您快速查閱！

CAT. NO.	DESCRIPTION	°C MIN/MAX	USE*
977200	OV-1 1% Chromosorb W-HP 80/100	100/350	General
977201	OV-1 3% Chromosorb W-HP 80/100	100/350	General
977202	OV-1 3% Gas Chrom Q 80/100	100/350	General
977203	OV-1 3% Chromosorb W-HP 100/120	100/350	General
977204	OV-1 3% Gas Chrom Q 100/120	100/350	General
977205	OV-17 3% Chromosorb W-HP 80/100	20/350	
977206	OV-17 3% Gas Chrom Q 80/100	20/350	General
977207	OV-17 3% Gas Chrom Q 100/120	20/350	Hydrocarbons
977208	OV-17 3% Chromosorb W-AW DMCS 100/120	20/350	Drugs Analysis
977209	OV-17 10% Chromosorb W-HP 100/120	20/350	Pesticides
977210	OV-25 3% Chromosorb W-HP 80/100	20/300	
977211	OV-101 1% Chromosorb W-HP 80/100	20/350	General
977212	OV-101 3% Chromosorb W-HP 80/100	20/350	Hydrocarbons
977213	OV-101 3% Gas Chrom Q 80/100	20/350	Carbohydrated
977214	OV-101 3% Chromosorb W-HP 100/120	20/350	FAME
977215	OV-101 10% Chromosorb W-HP 80/100	20/350	Alcohols
977216	OV-101 10% Chromosorb W-HP 100/120	20/350	
977217	OV-210 3% Chromosorb W-HP 80/100	20/275	General
977218	OV-210 10% Chromosorb W-HP 80/100	20/275	General
977219	OV-225 3% Chromosorb W-HP 80/100	20/250	General
977220	OV-225 3% Gas Chrom Q 80/100	20/250	General
977221	OV-225 10% Chromosorb W-HP 80/100	20/250	General
977222	OV-275 1% Chromosorb W-AW 80/100	20/250	
977223	OV-275 3% Chromosorb W-AW 80/100	20/250	General
977224	OV-275 3% Chromosorb W-AW 100/120	20/250	General
977225	OV-351 3% Chromosorb W-HP 80/100	50/250	General
977226	OV-351 10% Chromosorb W-HP 80/100	50/250	General
977227	Carbowax 20 M 3% Chromosorb W-HP 80/100	60/250	General
977228	Carbowax 20 M 10% + KOH 2% Chromosorb W-AW 80/100	60/230	Amines
977229	Carbowax 20 M 10% Chromosorb W-AW 80/100	60/250	General
977230	Carbowax 20 M-TPA 10% Chromosorb W-AW 80/100	60/250	Free Fatty Acids
977231	QF-1 3% Gas Chrom Q 80/100		General
977232	SE-30 GC Grade 3% Chromosorb W-HP 80/100	50/300	General
977233	SE-30 GC Grade 3% Gas Chrom Q 80/100	50/300	General
977234	SE-30 GC Grade 10% Chromosorb W-HP 80/100	50/300	General
977235	XE-60 3% Chromosorb W-HP 80/100	20/250	--
977236	Apiezon L 10% Chromosorb W-HP 80/100	20/100	--
977237	Dexsil 300 3% Chromosorb W-HP 100/120	20/450	HIGH TEMPERATURE
977238	DEGS 10% Chromosorb W-AW 80/100	20/200	FAME
977239	Tenax 60/80	0/375	--
977240	Tenax 80/100	0/375	--
977241	SE-30 4% + OV-210 6% Chromosorb W-HP 100/120	100/250	Pesticides
977242	OV-17 1.5% + QF-1 1.95% Gas Chrom Q 80/100	0/250	Pesticides
977243	OV-17 1.5% + QF-1 1.95% Chromosorb W-HP 100/120	0/250	Pesticides

★ Intended as a general guide.