



















# Injection Port Liners for Agilent Instruments

## Liners for Agilent Instruments

	Injection Type	ID x OD (mm)	Length (mm)	Packing	Cat. No. 5 Pack	Cat No. 25 Pack
	Direct Straight Liner	1.2 x 6.3	78.5	No	<b>453A1335</b>	<b>453A2332</b>
	Split Straight Liner	4 x 6.3	78.5	Quartz Wool	<b>453A2265</b>	<b>453A1262</b>
	Split Straight Liner	4 x 6.3	78.5	No	<b>453A1295</b>	<b>453A2292</b>
	Split/Splitless FocusLiner	4 x 6.3	78.5	Quartz Wool	<b>453A1255</b>	<b>453A1252</b>
	Split/Splitless FocusLiner with Single Taper	4 x 6.3	78.5	Quartz Wool	<b>453A1315</b>	<b>453A1312</b>
	Split/Splitless Liner with Single Taper	4 x 6.3	78.5	No	<b>453A1345</b>	<b>453A2342</b>
	Split/Splitless Liner with Double Taper	4 x 6.3	78.5	No	<b>453A1355</b>	<b>453A2352</b>
	Split/Splitless Liner w/ Recessed Gooseneck	4 x 6.3	78.5	Quartz Wool	<b>453A1305</b>	<b>N/A</b>
	Split/Splitless FAST FocusLiner	2.3 x 6.3	78.5	Quartz Wool	<b>453A1285</b>	<b>453A2282</b>
	Split/Splitless FAST FocusLiner with Single Taper	2.3 x 6.3	78.5	Quartz Wool	<b>453A2375</b>	<b>453A1372</b>
	Splitless Liner with Single Taper	4 x 6.3	78.5	Quartz Wool	<b>453A1925</b>	<b>N/A</b>
	Splitless Straight Liner	2.0 x 6.3	78.5	No	<b>453A2275</b>	<b>N/A</b>
	Splitless Liner with Recessed Gooseneck	2.0 x 6.3	78.5	No	<b>453A2325</b>	<b>N/A</b>
	Cyclo/Single Gooseneck (Deactivated Metal)	5.2 x 6.3	78.5	No	<b>453A2000</b>	<b>N/A</b>
	Single Gooseneck (Deactivated Metal)	5.2 x 6.3	78.5	No	<b>453A2001</b>	<b>N/A</b>
	Cyclosplitter Liner (Deactivated Metal)	5.2 x 6.3	78.5	No	<b>453A2002</b>	<b>N/A</b>
	Split/Splitless Liner with Wool (Deactivated Metal)	5.2 x 6.3	78.5	No	<b>453A2003</b>	<b>N/A</b>
	Split/Splitless Mixed Liner Sample Pack	Mixed	Mixed	No	<b>453AG001</b>	<b>N/A</b>

## Liner Sealing Rings for Agilent Instrument

Description	Cat. No.	Quantity
Graphite Liner Sealing Ring for Agilent Direct Straight Liners	<b>290GA242</b>	10 Pack
Graphite Liner Sealing Ring for Agilent SSL Liners	<b>290GA243</b>	10 Pack
Graphite Liner Sealing Ring for Agilent 2.0mm Liners	<b>290GA244</b>	10 Pack
Viton Liner Sealing Ring for Agilent SSL Injector	<b>290OA241</b>	10 Pack

Thermo Scientific GC consumables are suitable for use with Agilent GC systems

